





# Siron LIS (aPTT liquid)

REF	5035118	Siron LIS (aPTT liquid)	2 x 4 mL
REF	5035119	Siron LIS (aPTT liquid)	10 x 4 mL
REF	5035121	Siron LIS (aPTT liquid)	10 x 10 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola			
	manufacturer / Hersteller / fabbricante / fabricante / fabricant / Tilverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / Угробce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace / срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fídit se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определенний / роčet stanovení / Definicija
	distilled water / destilliertes Wasser / acqua distillata / agua destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / дестилирана вода / destilovaná voda / дестилированная вода / destilovaná voda / Serija		lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / партиа / партида номер / šarže / lot / šarže / in vitro dijagnostika
	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer		microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитърна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator		catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Конъюгат / Konjugát / Konjugat		ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk / έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola		stop solution / Stoppløsning / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stoppløsning / Stop-opløsning / Stoppløsning / διάλυμα πάυσης / Стоп разтвор / Стоп-разтвор / Zastavovací roztok / Stop solucija
	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o disolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplösa i / fortyndes eller oppløses i / Fortyndes eller oppløses i / αραιωση ή διάλυση σε / растворите или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / nafedte nebo rozpustte v / razrediti ili rastvoriti u		substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert / Inkubationsbuffer / Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer		washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solução de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösingskonsentrat / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миещ разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncetrat solucije za ispiranje
	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico en vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рηση διαγνωστικη εντόβ σωληνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnosticky prostfedek in vitro / Destilisana Voda		

**PRODUCT DESCRIPTION****INTENDED USE**

Siron LIS (lupus insensitive) is a liquid activated partial thromboplastin time (aPTT) reagent with **low sensitivity to lupus anticoagulants** to be used:

- as a screening test for the intrinsic coagulation system
- for specific determinations of factors VIII, IX, XI and XII when used in combination with the corresponding deficiency plasma.
- for control of heparin treatment.

**COMPOSITION**

Siron LIS is a liquid preparation of a stable aqueous suspension of phospholipids. FXII activation of Siron LIS is achieved by addition ellagic acid.

**MATERIAL REQUIRED** (not supplied with the kit)

- Buffer*			
<b>REF</b>	5277015	CaCl <sub>2</sub> 25 mmol/L	100 mL
- Control Plasma Normal and Abnormal *			
<b>REF</b>	5020040	Coagulation Control N	5 x 1 mL
<b>REF</b>	5021055	Coagulation Control A	5 x 1 mL

\*or any other package sizes, special Ceveron® alpha or TECHNOCLOT® Control and Calibration reagents of Technoclone.

**WARNING AND PRECAUTIONS**

- IVD for *in vitro* diagnostic use
- Siron LIS contains no human material. FXII activation of Siron is achieved by ellagic acid and phospholipids from natural origins.
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C.

Stability **opened** vials:

2...8°C 4 weeks	12°C Ceveron® alpha 10 days	18...25 °C 24 hours
--------------------	-----------------------------------	------------------------

**TEST PROCEDURE****PREPARATION OF PLASMA SAMPLES**

Plasma separation: Mix 9 parts of venous blood and 1 part of Sodium Citrate Solution (0.11 mol/L or 0.13 mol/L) and centrifuge for 15 minutes at a RCF of at least 2500 g (corresponding to DIN 58905).

**STORAGE OF PLASMA SAMPLES**

Plasma samples are stored in plastic tubes or siliconized vials at room temperature before testing:

- For patients that are not on heparin therapy, testing required within 3 hours
- For patients on heparin therapy, testing required within 1 hour
- For long term storage - 20°C, 2 weeks  
- 70°C, 6 months

**PREPARATION OF REAGENT**

Each vial contains 4mL / 10mL ready to use aPTT reagent. Before each use, let the reagent stabilize for 30 min at room temperature (18...25°C).

Some settling of the activator may occur during storage; thus, repeatedly invert Siron LIS gently just before use to ensure homogeneity. Continuous **stirring** is **not required** during use.

**PERFORMANCE OF THE TEST****CEVERON**

Technoclone provides application sheets for Ceveron® alpha. The application sheets contain analyser/assay specific handling and performance information which may differ from that provided in this instruction for use. In this case the information contained in the application sheets supersedes the information in this instruction for use. Please consult the instruction manual of the Ceveron® alpha.

**MANUAL**

The CaCl<sub>2</sub> 25 mmol/L should be preheated to 37°C.

Pipetting Scheme:

100 µL plasma sample
+ 100 µL <b>Siron LIS</b> in a prewarmed tube
incubate for <b>3 minutes</b>
+ 100 µL CaCl <sub>2</sub> solution (25 mmol/L; 37°C) Start stop-watch simultaneously
Leave the tube for 20 sec. before tilting. Determine the clotting end point when clot is formed. The results are given in seconds.

**ANALYSES RESULTS****REFERENCE RANGE**

**Normal range:** 30 - 44 seconds

The normal range refers to manual method. It is recommended that each laboratory establish and controls its normal range.

**Siron LIS in control of heparin therapy:**

A Siron LIS time of 2-4 times of the normal has been found to correspond to the accepted therapeutic level of heparin. However, the actual therapeutic level should be decided from the clinical situation.

**STANDARDISATION**

No international calibrator is available for the standardization of aPTT tests. Siron LIS documentation is based on studies testing plasma samples from normal subjects and heparinised plasma samples. An in-house reference batch is established in order to avoid batch-to batch variation.

**LIMITATION OF THE TEST**

- The blood samples should be taken in non-activating samples tubes in order to avoid activation of the clotting factors in the plasma sample.
- Plasma samples should be stored at room temperature to avoid cold activation.
- Platelet poor plasma samples should be used.
- It is important to use the correct incubation time, which is 4 minutes for Siron LIS when testing both heparinised and normal plasma samples.
- The calcium concentration used for recalcification should be 25 mmol/L when testing normal, abnormal or heparin plasmas.
- The reference range represents approximate values only. Each laboratory should determine its own normal range to take account of laboratory specific factors, such as the type of coagulometer used, reagent temperature, patient population, etc...
- The values found when testing control plasma should be compared to the values given on the control batch table for the corresponding lot of Control plasma. If the results obtained are outside the recommended range, avoid measuring patient samples until the problem is solved.
- A new calibration is required for each batch of Siron LIS and the additional reagents (e.g. deficient plasmas) where a calibration curve is necessary. In addition a new calibration is recommended, if software changes are introduced or after a major service of either instruments or equipment.
- Clotting time is not affected in this assay by samples containing ≤ 1000 mg/dL Haemoglobin, ≤ 46 mg/dL unconjugated bilirubin and ≤ 1000 mg/ml triglycerides. However elevated levels of conjugated bilirubin causes a prolonged clotting time.

**PERFORMANCE CHARACTERISTICS**

Performance data are given below. Results obtained in individual laboratories may differ.

**PRECISION**

Reproducibility was determined with different samples (in series and day to day). The following results were obtained:

Sample	Intra assay		Inter assay	
	Sample 1	Sample 2	Sample 3	Sample 4
n	21	21	15	15
MV (sec)	33.3	80.9	32.3	79.9
CV (%)	1.54	0.44	1.85	0.45

**COMPARISON OF METHODS OR CORRELATION**

Following correlation was obtained in comparing SIRON LIS with:

Actin® FS  $y = 1.04x - 2.77$   $R^2 = 0.986$

**APPLICATIONS FOR INSTRUMENTS**

Application Sheets are available from Technoclone or your local distributor upon request.

**LITERATURE**

Please contact Technoclone or your local distributor.

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