

Safety Data Sheet

Technoview Apixaban CON H



SDS No. 3D00.01
p 1/6 EN

Temp.: TF0366.01
Att 01 to D00068

ACCORDING TO EC-REGULATIONS 1907/2006 (REACH), 1272/2008 (CLP) & 2020/878

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product identifier

Product Name/Reference No. Technoview Apixaban CON H 5 x 1 mL REF 5090270

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified Use(s) For professional users only
Uses Advised Against Not known.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Company Identification Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH
Address of Manufacturer Brunner Strasse 67
1230 Vienna
Austria
Postal code 1230
Telephone: +43 1 86373-0
Fax +43 1 86373-44
E-mail products@technoclone.com
Office hours 8:00 – 16:00

1.4 Emergency telephone number

Emergency Phone No. +43 1 86373-0
Contact products@technoclone.com
National response center
Address NHS Direct
Emergency Phone No. +44 111

SECTION 2: HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture

Regulation (EC) No. 1272/2008 (CLP) Not classified as dangerous for supply/use

2.2 Label elements

Product Name According to Regulation (EC) No. 1272/2008 (CLP)
Technoview Apixaban CON H
Hazard Pictogram(s) None.

2.3 Other hazards

None known.

2.4 Additional Information

Contains human derived material - potentially biohazardous

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substances

Not applicable.

3.2 Mixtures

Product does not contain any hazardous ingredients.

SECTION 4: FIRST AID MEASURES

4.1 Description of first aid measures

Inhalation If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.
Skin Contact Wash skin with water.
Eye Contact Flush eyes with water for at least 15 minutes.

Safety Data Sheet



SDS No. 3D00.01
p 2/6 EN

Temp.: TF0366.01
Att 01 to DC0068

Technoview Apixaban CON H

Ingestion Wash out mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

None anticipated. Treat symptomatically.

4.3 Indication of any immediate medical attention and special treatment needed

Unlikely to be required but if necessary treat symptomatically.

SECTION 5: FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable Extinguishing media As appropriate for surrounding fire.
Unsuitable extinguishing media None.

5.2 Special hazards arising from the substance or mixture

None anticipated. Heating may cause decomposition.

5.3 Advice for firefighters

As appropriate for surrounding fire.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Wear suitable gloves if prolonged skin contact is likely.

6.2 Environmental precautions

Do not release large quantities into the surface water or into drains.

6.3 Methods and material for containment and cleaning up

Sweep spilled substances into containers if appropriate moisten first to prevent dusting.

6.4 Reference to other sections

See also Section 8, 13.

SECTION 7: HANDLING AND STORAGE

7.1 Precautions for safe handling

Not known.

7.2 Conditions for safe storage, including any incompatibilities

Storage temperature +2°C to +8°C
Storage life Stable under normal conditions.
Incompatible materials None known.

7.3 Specific end use(s)

In vitro diagnostic use

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

8.1.1 Occupational Exposure Limits No Occupational Exposure Limit assigned.

8.2 Exposure controls

8.2.1. Appropriate engineering controls Ensure adequate ventilation.

8.2.2. Personal protection equipment

Eye Protection Wear eye protection with side protection (EN166).
Skin protection Gloves recommended
Respiratory protection Normally no personal respiratory protection is necessary.

Safety Data Sheet



Temp.: TFF0366.01
Att 01 to D00068

Technoview Apixaban CON H

SDS No. 3D00.01
p 3/6 EN

Thermal hazards None known.
8.2.3. Environmental Exposure Controls Do not release large quantities into the surface water or into drains.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

Physical state	solid (lyoph.)
Colour	white-yellowish
Odour	odourless
Melting point/freezing point	Not known.
Boiling point or initial boiling point and boiling range	Not applicable.
Flammability	Not known.
Lower and upper explosion limit	Not known.
Flash Point	Not applicable.
Auto-ignition temperature	Not known.
Decomposition Temperature	Not known.
pH	Not known.
Kinematic Viscosity	Not known.
Solubility	Solubility (Water) : readily miscible. Solubility (Other) : Not known
Partition coefficient n-octanol/water (log value)	Not known.
Vapour pressure	Not known.
Density and/or relative density	Not known.
Relative vapour density	Not known.
Particle characteristics	Not known.

9.2 Other information

None.

SECTION 10: STABILITY AND REACTIVITY

10.1 Reactivity

None anticipated.

10.2 Chemical Stability

Stable under normal conditions.

10.3 Possibility of hazardous reactions

No hazardous reactions known if used for its intended purpose.

10.4 Conditions to avoid

None anticipated.

10.5 Incompatible materials

Not known.

10.6 Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity - Ingestion	Low acute toxicity
Acute toxicity - Skin Contact	Low acute toxicity
Acute toxicity - Inhalation	Low acute toxicity
Skin corrosion/irritation	Non-irritant
Serious eye damage/irritation	Non-irritant
Skin sensitization data	Not classified

Safety Data Sheet



Temp.: TFF0366.01
Att 01 to DC0068

Technoview Apixaban CON H

SDS No. 3D00.01
p 4/6 EN

Respiratory sensitization data	Not classified
Germ cell mutagenicity	There is no evidence of mutagenic potential.
Carcinogenicity	No evidence of carcinogenicity.
Reproductive toxicity	Not classified
Lactation	Not classified
STOT - single exposure	Not classified
STOT - repeated exposure	Not classified
Aspiration hazard	Not classified

11.2 Information on other hazards

Not known.

SECTION 12: ECOLOGICAL INFORMATION

12.1 Toxicity

Toxicity - Aquatic invertebrates	Low toxicity to invertebrates.
Toxicity - Fish	Low toxicity to fish.
Toxicity - Algae	Low toxicity to algae.
Toxicity - Sediment Compartment	Not classified.
Toxicity - Terrestrial Compartment	Not classified.

12.2 Persistence and degradability

Not known.

12.3 Bioaccumulative potential

Not known.

12.4 Mobility in soil

Not known.

12.5 Results of PBT and vPvB assessment

Not known.

12.6 Endocrine disrupting properties

None known.

12.7 Other adverse effects

Not known.

SECTION 13: DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Dispose at suitable refuse site.

13.2 Additional Information

No special precautions are required for this product.

SECTION 14: TRANSPORT INFORMATION

Not classified as hazardous for transport.

14.1 UN number or ID number

Not applicable

14.2 UN proper shipping name

Not applicable

14.3 Transport hazard class(es)

Not applicable

14.4 Packing group

Not applicable

Safety Data Sheet

Technoview Apixaban CON H



SDS No. 3D00.01
p 5/6 EN

Temp.: TF0366.01
Att 01 to DC0068

14.5 Environmental hazards

Not classified as a Marine Pollutant.

14.6 Special precautions for user

Not known

14.7 Maritime transport in bulk according to IMO instruments

Not known

SECTION 15: REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

European Regulations - Authorisations and/or Restrictions On Use

Candidate List of Substances of Very Not listed
High Concern for Authorisation

REACH: ANNEX XIV list of Not listed
substances subject to authorisation

REACH: Annex XVII Restrictions on Not listed
the manufacture, placing on the
market and use of certain dangerous
substances, mixtures and articles

Community Rolling Action Plan Not listed
(CoRAP)

Regulation (EC) N° 850/2004 of the Not listed
European Parliament and of the
Council on persistent organic
pollutants

Regulation (EC) N° 1005/2009 on Not listed
substances that deplete the ozone
layer

Regulation (EU) N° 649/2012 of the Not listed
European Parliament and of the
Council concerning the export and
import of hazardous chemicals

National regulations

Other Not known.

15.2 Chemical Safety Assessment

A REACH chemical safety assessment has not been carried out.

SECTION 16: OTHER INFORMATION

Date of Issue 2023-09-22 (replaces 356 from 2017-06-07)

Date of Revision n.a.

LEGEND

Hazard Pictogram(s) None.

Precautionary Statement(s) None.

Acronyms

ATE : Acute Toxicity Estimate
CAS : Chemical Abstracts Service
CLP : Regulation (EC) No 1272/2008 on classification, labelling and
packaging of substances and mixtures
DNEL : Derived No Effect Level
EC : European Community
EINECS : European Inventory of Existing Commercial Chemical
Substances
LTEL : Long term exposure limit
PBT : Persistent, Bioaccumulative and Toxic
PNEC : Predicted No Effect Concentration

Safety Data Sheet



Temp.: TFF0366.01
Att 01 to D00068

Technoview Apixaban CON H

SDS No. 3D00.01

p 6/6

EN

REACH : Registration, Evaluation, Authorisation and Restriction of Chemicals

STEL : Short term exposure limit

STOT : Specific Target Organ Toxicity

vPvB : very Persistent and very Bioaccumulative

Key literature references and sources for data used to compile the SDS

Regulation (EC) No. 1272/2008 (CLP)

Disclaimers

Information contained in this publication or as otherwise supplied to users is believed to be accurate and is given in good faith, but it is for the Users to satisfy themselves of the suitability of the product for their own particular purpose. Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition (statutory or otherwise) is excluded except to the extent that exclusion is prevented by law. Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH accepts no liability for loss or damage (other than that arising from death or personal injury caused by defective product, if proved), resulting from reliance on this information. Freedom under Patents, Copyright and Designs cannot be assumed.